## Transformation protocol

## Reaction 1 (cloning and transforming your gene)

- All reactions and tube must be on ice
- Prepare the ligation reaction mixture by combining (in order) the

following components:

1.

2.

6.

- 1.5  $\mu$  l StrataClone Cloning Buffer
- 0.5  $\mu$  l StrataClone Vector Mix amp/kan
- (Tube A containing 1.5  $\mu$  l Cloning and Buffer 0.5  $\mu$  l Vector Mix )
- 1  $\mu$  l of PCR product (gene purified from gel)
- 3. Mix gently by repeated pipetting, and then incubate the ligation reaction at room temperature for 5 minutes. When the incubation is complete, place the reaction on ice.
- 4. Add 1  $\mu$  l of the cloning reaction mixture to the tube of thawed competent cells. Mix gently (do not mix by repeated pipetting).
- 5. Incubate the transformation mixture on ice for 5–20 minutes.
  - Heat-shock the transformation mixture at 42°C for 45 seconds.
- 7. Incubate the transformation mixture on ice for 2 minutes.
- 8. Add 200  $\mu$  l of pre-warmed LB medium to the transformation reaction

mixture.

9. Allow the competent cells to recover for at least 1 hour at 37°C with agitation. (Lay the tube of cells on the shaker horizontally for better aeration.)

## Reaction 2 (cloning and transforming control gene)

- 1. Prepare the ligation reaction mixture by combining (in order) the following components.
- 1.5  $\mu$  l StrataClone Cloning Buffer
- 0.5  $\mu$  l StrataClone Vector Mix amp/kan
- 1  $\mu$  l of PCR product (gene purified from gel)
- (Tube B containing all three components)

- 2. Mix gently by repeated pipetting, and then incubate the ligation reaction at room temperature for 5 minutes. When the incubation is complete, place the reaction on ice.
- 3. Add 1  $\mu$  l of the cloning reaction mixture to the tube of thawed competent cells. Mix gently (do not mix by repeated pipetting).
- 4. Incubate the transformation mixture on ice for 5–20 minutes.
- 5. Heat-shock the transformation mixture at  $42^{\circ}$ C for 45 seconds.
- 6. Incubate the transformation mixture on ice for 2 minutes.
- 7. Add 200  $\mu$  l of pre-warmed LB medium to the transformation reaction mixture.
- 8. Allow the competent cells to recover for at least 1 hour at 37°C with agitation. (Lay the tube of cells on the shaker horizontally for better aeration.)

## Reaction 3 (transforming control vector)

(Tube C containing control puc18 vector)

- 1. Add 1  $\mu$  1 of the vector to the tube of thawed competent cells. Mix gently (do not mix by repeated pipetting).
- 2. Incubate the transformation mixture on ice for 5–20 minutes.
- 3. Heat-shock the transformation mixture at  $42^{\circ}$ C for 45 seconds.
- 4. Incubate the transformation mixture on ice for 2 minutes.
- 5. Add 200  $\mu$  l of pre-warmed LB medium to the transformation reaction mixture.
- 6. Allow the competent cells to recover for at least 1 hour at 37°C with agitation. (Lay the tube of cells on the shaker horizontally for better aeration.)