DNA extraction protocols

A. <u>DNA extraction from tissue using non-organic</u> <u>solvents:</u>

Buffers Preparation :

Sucrose buffer:

57 gm Sucrose

3.1 gm MgC12.6H2O

0.6 gm Tris. Hcl

500 ml Autoclaved d.d. H2O

adjust pH 7.5 with 0.1 N HCl. The solution is stable if stored at refrigerator for 1 week.

EDTA:

0.72 gm disodium EDTA250 ml Autoclaved d.d. H20adjust pH at 7.5 with 0.1 N Na0HStore at room temperatureEDTA inhibits the action of Dnase enzymes, makes thelysis of nuclear membrane easier.

SDS:

25 gm sodium dodecyl sulfate (SDS)
250 ml d.d. Autoclaved H2O
dissolve 2 gm SDS in 20 ml H2O
Store at room temperature
SDS emulsifies the plasma and nuclear membrane.

<u>2 M NaCl</u>

29.2 gm Nacl

250 ml Autoclaved d.d. H20

Store at room temperature.

Nacl increases ion concentration, which disrupt ionic bonds between DNA and Protein.

Protocol:

- 1) Dice a piece of fresh tissue "1 gm liver, 150 mg spleen".
- 2) Immerse diced tissue in isotonic sucrose buffer (10 ml for liver, 5 ml for spleen) in a blender and homogenize completely.
- 3) Strain homogenate through cheesecloth to remove clumps of gristle.
- 4) Pour homogenate equally in centrifuge tubes, spin for 10 min. at 8000 rpm.
- 5) Pour off suspension, resuspend the pellet in buffered sucrose (1.5 ml).
- 6) Transfer 1.5 ml well mixed suspended nuclei to small test tube.
- 7) Add 1 ml 28% EDTA, mix gently, let stand for 5 min.
- 8) Add 0.5 ml SDS, mix very gently and let stand for 2 min. Be careful at the following steps.
- 9) Add 300 µl prot. K incubate at 50°C for 30 min.
- 10)Slowly add 0.3 ml of 2 M Nacl Pause after each drop and mix gently.
- 11) Draw 2 ml of 95% chilled ethanol, eject it slowly down the side of test tube.
- 12) Using a long stirring rod, make a slow stirring motion and spool the protein free DNA.
- 13) DNA may be saved in ethanol resuspended in one ml TE buffer.

B. <u>Extraction of Chromosomal bacterial DNA by</u> <u>diatoms (Total or genomic bacterial DNA):</u> A) Preparation of diatoms suspension:

- 1) Add 50 ml H2O and 500 µl of 35% Hcl to 10 g celite or diatoms.
- 2) Divide the suspension into aliquots into small bottles which are tightly closed and autoclave, keep in room temperature in the dark.

B) <u>Buffers:</u> Lysis buffer (L6)

Dissolve 30 gm GuScN into 25 ml of 0.1 M Tris Hcl (pH 6.4) + 5.5 ml of 0.2 M EDTA. Adjusted with NaOH to pH 8.0 + 0.65 ml Triton X-100. Homogenize the solution in shaking water bath at 60°C (stable at room temp. in the dark for at least 3 weeks).

Washing buffer (L2)

Dissolve 30 gm of GuscN in 25 of 0.1 M Tris Hcl (pH 6.4) by heating in shaking water bath at 60 - 65°C. (stable at room temp. in the dark for at least 3 weeks).

Elution buffer (TE)

- 10 mM Tris-Hcl
- 1 mM EDTA pH 8
- 100 X concentrated stock solution must diluted to 1 X working solution.
- These buffers prepared under fume hood.

Reaction vessels

Add 900 μ l of lysis buffer (L6) and 40 μ l of diatoms suspension (shake well before use) to 1.5 ml ep-tube homogenize the solution by vortexing.

Reaction vessels can be successfully used up to 1 week after assembly.

C) Method:

- 1) Vortex the reaction vessel.
- 2) Vortex the overnight bacterial broth culture and take 1 ml in eppindorf tube or 50µl sample (urine or serum).
- 3) Centrifuge at 13000 for 15 min, then bacterial pellet is suspended in 50 μ l autoclaved distilled water or TE buffer, vortex.
- 4) Add the suspended pellet to the reaction vessel.
- 5) Leave at room temperature for 10 min. vortex several times during the incubation time.
- 6) Centrifuge for 2 min. at 13000 rpm.
- 7) Dispose the supernatant.
- 8) Diatom NA pellet is washed twice with washing buffer L2 then twice. with ethanol 70% and once with acetone. Vortex then centrifuge for one minute each time.
- 9) Dry the vessel at 56°C with open lid ep-tube (use the thermal cycler **or** eppendorf heat block) for 10 min.
- 10)Add elution buffer TE with or without Rnase then vortex briefly and incubate for 10 min at 56°C (100-75 μ l TE).
- 11)Vortex again and centrifuge at 13000 for 3 min.
- 12)Use the supernatant containing DNA for further experiment.

C. Extraction of Plasmid DNA (Alkaline lysis		
<u>method):</u>		
<u>Preparation of media & solutions</u>		
* <u>L B Broth (</u> pH_7.0)		
- Tryptone	10 gm	
- Yeast extracte	5 gm	
- Nacl	5 gm	
- 4 N NaoH	1-2 ml	
- D.W.	1000 ml	
* <u>Sol.1</u>		
- 50% glucose solution	1.8 ml	
- 250 mM EDTA - pH8	4.0 ml	
- 1 M tris pH 8.1	2.5 ml	
- D.W.	60 ml	
* Sol II		

a) 10% SDS	store	-20 in aliquots
b) 2 N NaoH	store	-20 in aliquots

Just before use, mix 1 ml of (a) + ml of (b) + 8 ml D. Water (the total volume is 10 ml).

Electrophoresis

- 0.7% agarose
- 10 ml of DNA + 5 μ l dye
- 1µl marker + 7µl TE+ 5µl dye

<u>Method:</u>

1) Inoculate single colony in 5 ml LB Broth.

2) Incubate at 37°C for 10 hour (over night) in shaking bath.

- 3) Take 2 ml of the culture in eppendorf tube centrifuge for 1 min at 15,000 rpm at 4°C.
- 4) Aspirate supernatant and take the pellet.

- 5) Add 100 μ l of solution (I) and vigorous vortexing till the pellet disappear.
- 6) Add 200 μ l of solution 11 without vortex or shaking, just gentle inversions. keep on ice bath for 5 min.
- 7) Add 150 µl of 3 M sodium acetate (pH 4.8), gentle inversion and incubate on ice for 5 min.
- 8) Centrifuge 5 min at 15000 rpm at 4°C.
- 9) Transfer clear lysate to ep-tube.
- 10)Add double volume of cold ethanol with several inversions, incubate at-20 for 10 min.
- 11) Centrifuge for 5 min at 15,000 at 4°C.
- 12) Dry the pellet for 5-10 min.
- 13)Resuspend in 40 μ l of TE pH 8.

14) Store at -20.

D- DNA extraction from blood:-

Procedure:

- 1) Add 1 ml of the DNA extraction buffer (provided with the kit), mix well till complete dissolve of the pellet and incubate 30 min in water bath at 65°C.
- 2) Spin down for 15 min (13000 xg) and transfer the upper 500 μ l of the supernatant to a new tube containing 500 μ l ice-cold isopropanol.
- 3) Mix well and wait 5-25 min at RT (you can wait longer).
- 4) Spin down for 20 min (13000 xg).
- 5) Gently, pour off the supernatant.
- 6) Wash pellets two times with 500 μ l 70% EtOH for 5 min and air dry.
- 7) According to the size of the resulted pellet, resuspend pellets in 50-100 µl preheated (55°C) sterile water or 1x TE buffer {10 mM Tris, 1 mM Na2 EDTA (pH 8)} and then load 5-10 µl on Agarose gel to check.
- 8) Store at -20°C until use.