An Introduction to Polymerase Chain Reaction (PCR)




Prof. Dr. Hamdy M. El-Aref Assiut University, Genetics Department, Faulty of Agriculture



Introduction

- > The technique was invented by Dr. Kary Mullis, 1983,
- For which he received the Nobel Prize in Chemistry in 1993.



in vitro method (test tube system) for DNA replication.



Published papers with 'PCR'

1989 - 219 1990 - 4961991 - 711 1992 - 9061993 - 10301994 - 857 (>4000) 1995 - 8231996 - 7961997 - 7322006,3,22 - 255,788 2006/4/18 - 257,737 2007/3/9 - 283,607 2007/4/11 - 286,486

1998,10 - >73,000 1999,4 - >81,000 2000,10 - 121,305 2001,2 - 125,563 2002,3 - 149,572 2003,2 - 170,841 2004,2,23-195,193 2004,2,26-195,265




Introduction

PCR will allow a short stretch (specific sequence) of DNA (usually fewer than 3000 base pairs) to be amplified to more than million fold.

Millions of copies of a segment of DNA can be made within a few hours.

(n) nitro^{PDF} professiona




download the free trial online at nitropdf.com/professional

What do we need for PCR?

- Target DNA (Template).
- Two primers: (forward and reverse)
- Nucleotides: (the 4 dNTP'S: A, T, C, G)
- Heat-stable DNA polymerase:
 - (like Taq DNA Polymerase)
- Buffer and Cofactor MgCl₂ (Mg⁺⁺, K⁺)⁻
 - Thermal cycler.





The Reaction

PCR Eppendorf tube containing reaction mix. placed in a thermocycler.







PCR tube

THERMOCYCLER

Created with

nitro^{PDF} professional


PCR Cycle


Each cycle (Round) of PCR contains 3 steps:

1- Denaturation


lab

- 2- Primer annealing
- **3- Primer extension**
- The cycle usually repeated for 25 40 times.







download the free trial online at nitropdf.com/professional




download the free trial online at nitropdf.com/professional






The mixture is heated again at 90-95°C to denature the molecules and separate the strands and the cycle is repeated.

Each new strand then acts as a template for the next cycle of synthesis.





download the free trial online at nitropdf.com/professional

A thermal cycler automatically changes the temperature at the correct time for each of the stages.

How much amplification do you get? • Amplicons = A X 2^{n-2}

n = number of cycles

>A = starting target copy number


Ex.
 6 cycles = 26⁶⁻² (or 16) DNA copies.
 26 cycles = 2²⁴ (or 16.77 million) DNA copies.

PCR animation and Video



The amplified product can be detected using gel electrophoresis to view the band containing DNA fragments.







Primer features:

single stranded DNA molecules.

•Types of primers

- Random (short 10 15 b) [e.g. RAPD-PCR]
- Specific (long 17 24 b)
 - *as the primer increases in size, the chances of matching the target size increase.

*Longer the primer, the higher the annealing temperature.

They are manufactured commercially

Tag DNA Polymerase features :

A heat-stable DNA polymerase must be used in the reaction.
Heat stable Tag DNA polymerase

Heat stable <u>Taq DNA polymerase</u> was purified from the hot springs (hydrothermal vents) bacterium *Thermus aquaticus* in 1976.

Taq polymerase has optimal enzymatic activity at 72°C.



Created with



Tag DNA Polymerase features :

Its enzymatic halflife (at 95°C) is 40 min

Taq Polymerase extends the DNA chain by adding ~1.0 Kb per min.




Polymerase	T ½, 95℃	Extension Rate (nt/sec)	Source
Amplitaq	80 min	>50	T. aquaticus
Vent	400 min	>80	Thermococcus litoralis
Deep Vent	1380 min	?	Pyrococcus GB - D
Pfu	>120 min	60	Pyrococcus furiosus
Tth	20 min	>33	T.

Thermal Cycler:

It heats and cools the reaction tubes within it to the precise temperature required for each step of the reaction.

It is available in different specificity and models.






download the free trial online at nitropdf.com/professional

Programming the Thermocycler

otect Lab



download the free trial online at nitropdf.com/professional

าal

Advantages of PCR

- PCR is fast (2–5 hours).
- DNA or RNA can be amplified.
- High-yield amplification can be achieved (10⁶ to 10⁹ amplification).
- DNA from one cell equivalent can be amplified.

• Works on damaged DNA.

- PCR products can be directly sequenced.
- Flexible.



Disadvantages of PCR

- Need information about Target DNA sequence.
- Highly susceptible to contamination or false amplification.
- Amplification may not be 100% specific.
- Specificity of amplification is dependent on temperature and Mg++ concentration.






Disadvantages of PCR

 Analysis and product detection usually takes longer than the PCR reaction itself.

There is an upper limit to the size of DNA synthesized by PCR.





References

- Micklos, David, Greg Freyer and David Crotty. DNA Science a First Course. New York:Cold Spring Habor Laboratory Press, 2003.
- Purves, Sadava, Orians, Heller. "Life." 6th ed. Sinauer Associates, 2001.
- Demidov.V, Broude. N(2004). DNA Amplification: Current Technologies and Applications

• Websites:

- http://faculty.plattsburgh.edu/donald.slish/PCRmov.html (Animation)
- http://www.accessexcellence.org/RC/AB/IE/PCR_Xeroxing_DNA.html
- http://www.people.virginia.edu/~rjh9u/pcranim.html (PCR Animation)
- http://www.escience.ws/b572/L3/L3.htm (PCR Animation)
- http://homepages.strath.ac.uk/~dfs99109/BB211/RecombDNAtechlect4.html
- http://en.wikipedia.org/wiki/Polymerase_chain_reaction
- http://www.escience.ws/b572/L3/L3.htm
- http://allserv.rug.ac.be/~avierstr/principles/pcrani.html

nitro^{PDF} professional

